



Which wireless solution is right for your patient?



2.4 GHz TECHNOLOGY

900SYNC™ TECHNOLOGY


PRODUCTS & STYLES



WIRELESS



CONNECTIVITY



KEY FEATURES

2.4 GHz & 900SYNC TECHNOLOGY

- > Industry-best feedback cancellation
- > Best-in-class adaptive noise management
- > Ear-to-ear synchronized user controls
- > Continuously adaptive directionality
- > Smart frequency-lowering technology
- > Multiflex Tinnitus Technology
- > Dedicated music compression focused on listening enjoyment
- > Wireless CROS/BiCROS¹





2.4 GHZ ACCESSORY GUIDE

		APPLE	ANDROID	MINI REMOTE MICROPHONE	TV	REMOTE MICROPHONE +	REMOTE
Bluetooth® Technology	Mobile Phone Call Streaming	●	● ²			●	
	Media Streaming	●	● ²			●	
SNR Enhancement and Comfort Features	Surround Remote Microphone					●	
	Focus Remote Microphone			●		●	
	On-Demand Noise Control	●	●				●
	User-Created Custom Memories	●	●				
Personalization Tools	Automatic Memory Engagement (Geotag)	●	●				
	Noise Control	●	●				
	Automatic Car Memory	●	●				
	Automatic Phone Memory	●	●				
	Multiflex Tinnitus Adjustments	●	●				
Audio Streaming	TV Streaming			●	●	●	
	Line-In Audio Streaming from TV, Computer, MP3, etc.				●	●	
Hearing Aid Controls	Hearing Aid Volume Changes	●	●				●
	Hearing Aid Memory Changes	●	●				●
	Mute	●	●				●
	Home Shortcut						●
Convenience Features	Multiflex Tinnitus ON/OFF	●	●				●
	Find My Hearing Aids	●	●				
	Accessory Volume	●	●		●	●	●
	Accessory Start/Stop	●	●				●
	Teleloop					●	
	FM Receiver/DAI					●	

● Standard feature ● Feature available with accessory indicated



900SYNC ACCESSORY GUIDE

		MINI MOBILE	REMOTE MICROPHONE 2	MEDIA DEVICE 2	REMOTE CONTROL	
					ADVANCED	BASIC
Bluetooth Technology (iOS and Android Compatible)	Mobile Phone Call Streaming (Apple and Android)	●				
	Audio Streaming (Apple and Android)	●				
SNR Enhancement and Comfort Features	Surround Microphone	●	●			
	Focus Microphone	●	●			
	On-Demand Noise Control				●	
	TV Streaming	●	●	●		
Audio Streaming	Line-In Audio Streaming from TV, Computer, etc.	●	●	●		
	Volume Changes			●	●	●
Hearing Aid Controls	Memory Changes				●	
	Mute				●	
	Home Shortcut				●	●
Convenience Features	Multiflex Tinnitus ON/OFF				●	

¹ Available on all standard products except the Livio and Livio AI micro RIC 312.

² Direct streaming with Google Pixel 3, 3XL, 3a, 3aXL. Remote Microphone + required for streaming with other Android devices.

For Livio AI or Livio, Thrive Hearing Control app and iOS/Android compatibility information, please visit starkey.com/thrive-hearing.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Starkey is under license.

Use of the Made for Apple badge means that an accessory has been designed to connect specifically to the Apple product(s) identified in the badge, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a registered service mark of Apple Inc.

Android and Google Play are trademarks of Google Inc.

Galaxy is a registered trademark of Samsung Electronics Co., Ltd.

iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.