

# Rechargeable Shipping Guide

If you are sending hearing aid(s) and charger(s):



This applies if hearing aid(s) and charger(s) are being sent with OR without 1 or more accessory.

## Labels Needed:

- UPS Overnight Label
- Batteries Contained Label – UN3481



Do not print

## Instructions:

- 1 Print labels to scale AND in **color**
- 2 Cut individually and securely adhere to box along with the shipping label
- 3 Labels must be on the same side of box and cannot touch or be folded over any side of the box
- 4 Items must be turned off and secured against movement, damage or accidental activation
- 5 Items must be secured in a box



## If you are sending charger(s) only:



This applies if 1 or more charger is being sent in without the hearing aids.

### Labels Needed:

- UPS Ground Label
- Batteries Contained Label – UN3480
- Lithium Ion Battery Label



Do not print

**LITHIUM ION BATTERIES -  
FORBIDDEN FOR  
TRANSPORT ABOARD  
PASSENGER AIRCRAFT**

9/20\_86478-007

Do not print

## Instructions:

- 1 Print labels to scale AND in **color**
- 2 Cut individually and securely adhere to box along with the shipping label
- 3 Labels must be on the same side of box and cannot touch or be folded over any side of the box
- 4 Chargers must be packed individually inside shipping box
- 5 Outer package must be a box

## If you are sending 1 or more accessory:



This applies if 1 or more accessory is being sent with OR without the hearing aids.

*Note: If there are no hearing aids in the shipment, make sure to send any charger in a separate package.*

### Labels Needed:

- UPS overnight label
- Batteries Contained Label – UN3481



Do not print

## Instructions:

- 1 Print labels to scale AND in **color**
- 2 Cut individually and securely adhere to box along with the shipping label
- 3 Labels must be on the same side of box and cannot touch or be folded over any side of the box
- 4 Items must be turned off and secured against movement, damage or accidental activation
- 5 Items must be secured in a box

## Printable Battery Labels

Label needs to be printed to scale. In print menu, select “Custom Scale: 100%”.



Batteries Contained Label – UN3481

## Printable Battery Labels

Label needs to be printed to scale. In print menu, select "Custom Scale: 100%".



Batteries Contained Label – UN3480

**LITHIUM ION BATTERIES -  
FORBIDDEN FOR  
TRANSPORT ABOARD  
PASSENGER AIRCRAFT**

9/20\_86478-007

Lithium Ion Battery Label